



Medical Disposables Corp.  
YES, WE HAVE IT. IT'S IN OUR NAME.

# MD Fentanyl

## Single Test DipCard

99% Accuracy & Fast Results in 5 minutes

### PRODUCT FEATURES:

- CLIA Waived / OTC
- Built-In controls
- Non-invasive collection
- 10-15 seconds dipcard
- Results in 5 minutes
- Ultra-Sensitive detection 1.0 ng/mL
- Results can be photocopied
- Accurate urine rapid test dipcard



Item Code: MDFN-114C

**Ultra-Sensitive Detection**

Contact Us at [Sales@MedicalDisposables.us](mailto:Sales@MedicalDisposables.us)

### MD CLIA Waived OTC Fentanyl Single Test Card

The MD CLIA Waived OTC Fentanyl Single Drug Test Card addresses the Opioid crisis by offering a fast, reliable, and accessible solution for detecting fentanyl in human urine. With *ultra-sensitive detection* at a low cutoff concentration of **1.0 ng/mL**, it ensures no trace goes unnoticed. This CLIA-waived and over-the-counter (OTC) product is approved for simple, non-prescription use in professional, medical, and home settings. Delivering accurate results in just 5 minutes and requiring no special equipment, this portable test card is the ideal tool for on-site testing and rapid decision-making.



(888) 863-1112

6831 Edgewater Commerce Pkwy Ste 1101  
Orlando, FL 32810

MedicalDisposables.us  
[Sales@MedicalDisposables.us](mailto:Sales@MedicalDisposables.us)

Like Us!



Follow Us!



Follow Us!




**HAVE QUESTIONS? PLEASE CALL US (888) 863-1112**

# Single Urine Test Dipcard

## Step 1

Collect urine from the donor in a disposable collection cup\*. If the collection cup has a temperature strip read it and make sure is between temperature range.

\* Sold separately




## Step 2

Remove the test device from the foil pouch. Remove the cap from the test device. Label the device with patient or control identifications.



## Step 3

Immerse the absorbent end into the urine sample for about 10 - 15 seconds. Make sure that the urine level is not above the marker line printed on the front of the device.



## Step 4

Lay the device flat on a clean, dry, non-absorbent surface.

## Step 5

Read the results in 5 minutes. Do not interpret the results before 5 minutes or after 10 minutes, as they may be invalid.



## INTERPRETING RESULTS

### NEGATIVE



Two lines appear. One red line should be in the control region (C), and another apparent red or pink line adjacent should be in the test region (T). This negative result indicates that the drug concentration is below the detectable level.

### POSITIVE



One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the drug concentration is above the detectable level.

### INVALID



Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact your Medical Disposables Corp. representative.



(888) 863-1112

6831 Edgewater Commerce Pkwy Ste 1101  
Orlando, FL 32810

MedicalDisposables.us  
Sales@MedicalDisposables.us

Like Us!



Follow Us!



Follow Us!



HAVE QUESTIONS? PLEASE CALL US (888) 863-1112