

Disposable
Surgical Cautery Pen



EndoGuidance[®]
ELECTROSURGICAL PENCIL

YES, WE HAVE IT. IT'S IN OUR **NAME**[®]

Surgical Division

More **Products**
More **Solutions**
Solid **Results**





Electrosurgical Disposable Cautery Pen

Transform surgical precision with multi-functional operational dissectors. The ultimate tool for efficiency, precision, and versatility in surgery. Designed for flexibility, these dissectors enable surgeons to switch between functions effortlessly, without compromising precision or control. This innovation minimizes surgical time and maximizes patient safety.

Versatile Functions:

Dissect, cut, and coagulate with a single tool, reducing the need for multiple instruments during procedures.

Advanced Cutting Technology:

Ensures clean and controlled cuts, reducing tissue trauma and promoting faster healing.

Sterile and Safe:

Designed with patient safety in mind, our dissectors are fully sterilizable and meet stringent medical standards.

Ergonomic Design:

Engineered for comfort, allowing surgeons to maintain precision even during lengthy surgeries.

Easy Integration:

Compatible with most standard operating room systems and setups, making adaptation and integration simple.

Latex Free:

Our latex-free Electrosurgical Pencils ensure safety by eliminating allergy risks, making them ideal for all surgical environments.

REF. MD-ESU-D006

1 Box (25 units)

1 Carton (50 units)

The design of the texture allows the doctor to have better maneuverability in various operations.

The Cut and Coag Buttons provide easy control for precise cutting and coagulation during procedures.

The 3-Pin Plug ensures secure, reliable connectivity with electrosurgical units for optimal performance.

The special process coating treatment (PTFE) of the tip reduces eschar and produces less smoke.



Disposable Electro Surgical Pencil
Single use to avoid cross-infection



Medical Disposables Corp.
YES, WE HAVE IT. IT'S IN OUR NAME

Like Us!



Follow Us!



Follow Us!



(888) 863 -1112

6831 Edgewater Commerce Pkwy Ste 1101
Orlando, FL 32810

MedicalDisposables.us
SurgicalDivision@MedicalDisposables.us

HAVE QUESTIONS? PLEASE CALL (888) 863-8112