

Disposable  
Sterile Disposable Skin Stapler




**EndoGuidance**<sup>®</sup>  
SKIN STAPLER

YES, WE HAVE IT. IT'S IN OUR **NAME**<sup>®</sup>

*Surgical Division*

*More **Products**  
More **Solutions**  
Solid **Results***





Laparoscopic

# OBS® Sterile Disposable Skin Stapler

The OBS® Sterile Disposable Skin Stapler ensures precision and efficiency in suturing, making it an essential tool across diverse surgical disciplines. With its user-friendly design and exceptional performance, it supports successful outcomes in even the most delicate surgical procedures.

### Ergonomic Design:

Easy and convenient to operate with a concise, intuitive design.

### Precision & Safety:

Ensures a smooth, efficient procedure with accurate needle withdrawal.

### Aesthetic Results:

Produces a smooth and beautiful suture with minimal scarring.

### Biocompatibility:

Excellent tissue compatibility with reduced inflammatory response.

### Quick Recovery:

Easy and convenient stitch removal promotes faster wound healing.


### Enhanced Patient Comfort:

Reduces pain and improves operational efficiency.


Product Code	OBS-MW
Size	W: 7.0mm x 4.0mm Pin diameter 0.6mm
10 Box + 1 OBS® Disposable Staple Remover	

It features a slim, lightweight body that fits most hand sizes. The unique angled head provides a clear view of the work site from your normal vantage point.

The unique angled head provides a clear view of the work site from your normal vantage point.



The clear side window lets you view the number of staples remaining to minimize wastage.



Each stapler is preloaded with 35 wide stainless-steel staples.



**Medical Disposables Corp.**  
YES, WE HAVE IT. IT'S IN OUR NAME

Like Us!



Follow Us!



Follow Us!



(888) 863 -1112

6831 Edgewater Commerce Pkwy Ste 1101  
Orlando, FL 32810

MedicalDisposables.us

SurgicalDivision@MedicalDisposables.us

**HAVE QUESTIONS? PLEASE CALL (888) 863-8112**